
UR-WS01

General guidelines

- Take the time to thoroughly read the following instruction and utilize the product accordingly.
- Ensure you retain this manual and pass it along when transferring the product.
- Please be aware that this summary may not encompass every detail of all variations and deliberated steps. Feel free to reach out to us for additional information and assistance.

Notes

- The maximum static load capacity is 100 kg (220 lbs), this must not be exceeded.
- Before assembly, carefully check the parts for damage. Improper assembly or damaged parts pose a risk of injury.
- Before using the walking stick, it is advisable to have a specialist, or GP, look at it and help adjust it to your height.
- Under no circumstances use the product if any part is missing or defective.
- Note: Do not use the walking stick on wet and slippery surfaces. You will risk slipping and injury.
- Do not hang bags or other objects on the walking stick while using the walking stick. This will negatively affect its stability during use.
- Regularly check the strength of the joints.
- Try to be as considerate as possible when using the walking stick on your route.
- Never use the walking stick for any purpose other than support.
- Take extra care when using the walking stick on steep sections, stairs, thresholds and pavements.
- This product is intended for indoor and outdoor use by one person.
- To set the correct height, the handle of the walking stick should be at the level of your wrist in a natural, relaxed position with the arms alongside the body.
- When unfolding the walking stick, it is very important to check the joints before using the walking stick and leaning on it to avoid injury.

Warnings

- Keep all parts and packaging materials (such as foil, plastic bags and foam) away from babies and children to avoid choking hazards.
- Keep small parts out of the reach of children during assembly. Swallowing or inhaling these parts can be fatal.
- Always use the walking stick properly. It is strictly forbidden to use the walking stick as a climbing aid.
- Children must be supervised by adults when using the product.
- Changes or modifications to the product without permission may compromise safety.

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